



# Etac Prio 3A

Customise your wheelchair  
with accessories and options

 **etac**<sup>®</sup>

# Accessories & options

Customise the wheelchair Etac Prio 3A with a wide selection of accessories and options to choose from.

Here we list a selection of possibilities. More accessories and options can be found on [etac.com](http://etac.com) and in configuration.



### Back cushion

Soft and flexible back cushion, attached with velcro to be combined with side cushions and wedges.

Cover: Grey plush or grey hygiene fabric



### Seat cushion, Basic

The Basic cushion is gently contoured at the rear.

Cover: Grey plush or grey hygiene fabric



### 3D back cushion

Unique back cushion providing support and space for a kyphotic back or asymmetries, attached with velcro. Can be combined with side cushions and wedges.

Cover: Grey plush



### Seat cushion, Comfort

The comfort cushion is slightly lower at the rear. To further improve the stability, it has a removable viscoelastic foam insert and contouring for the thighs.

Cover: Grey plush or grey hygiene fabric.



### Back support Cover

A thin cover that provides a wider back. Can be combined with side cushions and wedges.

Cover: Grey plush



### Solid seat

Available in wood or plastic.



### Dynamic posterior pelvic support

Reclining the back support will automatically reduce the strap tension and retighten when the back is brought to an upright position. The ideal combination for activity and resting. The pelvic strap is height adjustable for precise setting.



### Seat Velcro attachment kit

To be fastened on the solid seat and on the underside of the seat cushion.



### Back cushion extender

Extend the back cushion when a higher back is used to cover the gap between the seat and the back cushion and increase the support surface.

Cover: Grey hygiene fabric



### Head support

Adjustable in height, depth and sideways. The cover is detachable and washable.

Cover: Grey plush or grey hygiene fabric



### Extension back upholstery

Extension part for the back upholstery to discreetly cover the gap between the seat and the back at higher back heights, for example if the user uses incontinence protection. The back upholstery also provides extra support for the seat cushion.



### Head support with side support

Adjustable in height, depth, angle and sideways. The cover is detachable and washable.

Cover: Grey hygiene fabric



**R82, anatomic head support**

Anatomically shaped with side supports. Padded and with a comfortable cover. Adjustable in height, depth, angle and sideways. The cover is detachable and washable.



**Memory function**

Head support can be equipped with a memory function to ensure correct installation after removal and replaced.



**Arm support, std**

Height and depth-adjustable, detachable.

Length: 49 cm (19 3/4")  
Height: 23,5-35 cm (9 1/4-13 3/4")



**Arm support, narrow**

Provides better access to the wheels when self propelling. Height adjustable and detachable.

Length: 38 cm (14 11/12")  
Height: 21-32,5 cm (8 2/4-12 2/3")



**Arm support soft covers**

Soft padding is available for our arm supports, standard and narrow.

Cover: Grey plush



**Arm support pad, hemi**

Adjustable in angle, the height is set with the arm support adjustment. Fitted onto the arm support narrow.

Cover: Padded detachable in grey plush



**Upholstery sideguard**

Reduces the seat width with 2,5 cm (1"). To be attached on the side guards. Sold per piece.

Cover: Grey hygiene fabric



**Tray**

Transparent, mounted on the arm support.



**Hemi Tray**

Fastened to the arm support cushions using Velcro (included). Fits both short and narrow arm support. Available for seats up to 50 cm (20") width. A non-slip mat is available as an option.



**Widening Kit**


Increases the width between the arm supports with 2 cm (3/4") by adding spacers.



**Leg support, elevating**

Swing-away and detachable. With foldable calf supports, adjustable in height and depth.

Cover: Grey hygiene fabric



**Leg support, Cross**

Swing-away and detachable. Available in 90° or 70°. In combination with 90° we recommend wide front castor housing.



**Foot support**

Height and angle adjustable. Can be rotated for various depth installation.  
Length: 13 cm (5").



**Heel strap**

Adjustable length. Cannot be combined with the extended foot support.



**Add-On Extended foot supports one piece**

One platform that is mounted on to existing foot supports. Available in one size:  
Small - 250 x 225 mm  
*Use the widened castor housing to avoid a collision.*



**Calf strap padding**

Detachable. Available in lengths of 20, 30 and 40 cm.



**Add-On Extended foot supports**

Is mounted on to existing foot supports. Use separately or bolt together and use as platform. Available in two sizes:  
Medium 150 x 225 mm  
Large 200 x 225 mm  
*Use the widened castor housing to avoid a collision.*



**Padded Calf strap**

Soft and warm comfort, Velcro adjustment. Assembly on fixed or elevating leg support. Available in 2 sizes: 70 cm and 90 cm



**Add-On Side protection**

The side protection is mounted on to the Add-On Extended footplates, to prevent the foot from coming in to contact with the castor. Sold per piece.



**Add-On Distal thigh support**

Helps to align the femurs, position the knees in neutral alignment and to control abduction.



**Add-On Footbox**

Is mounted on to existing foot supports Available in two sizes:  
Medium, for seat width 40 - 45 cm  
Large, for seat width 50 - 52,5 cm  
*Use the widened castor housing to avoid a collision.*



**Leg abduction**

Detachable and depth-adjustable.  
Cover: Grey hygiene fabric



**Extended foot support**

Fit the Cross foot support.  
Length: 21 cm (8 3/4").



**Amputee support**

Adjustable in height, depth, angle and sideways. Fits both left and right side. Available in plush and hygiene fabric.



**Foot support cover**

For added comfort, plush



**Trunk support, swing-away**

Adjustable in height, depth and width. The arm support can be removed when the lateral is swung out. Ideal for side transfer.

Cushion: Polyurethane



**Calf strap**

Length-adjustable and detachable.



**Trunk support, adjustable**

Adjustable in height, depth and width. Easily removed and installed without tools.

Cushion: Polyurethane



### Cover trunk support

Cover for trunk support in grey plush with a small straight cushion.



### Side Cushion

Embracing support. Rotate to choose between two different shapes or to be combined with Side stop. Sold per piece.

Cover: Grey plush



### Side stop, straight, 3A

The side stop is mounted on to the 3A back support. Combine with wedges or side cushion. Adjustable in depth and height. Sold per piece.



### Pressure distributor

Distributes pressure from the back profile. Can also be used as a widening of the back. Sold per piece.



### Width extension back 30 mm, 3A

The Width extension is mounted on to the 3A back support, and extends the width by 30 mm (1") per side. Combine with wedges, pressure distributors and larger cover. Sold per piece.



### Back wedges

Soft wedge that attaches with velcro under the cover. Combine with side stop straight or with width extension. Sold per piece. Two sizes: 15 cm (6") and 30 cm (12") long. Cover: Dark grey velour.



### Wedges cell foam

Cell foam with velcro attachment. Sold in pair. Available in three sizes:  
Length x Width x Height  
Small: 13 x 9,5 x 2 cm (5 x 4 x 3/4")  
Medium: 19 x 12 x 2 cm (7 1/2 x 5 x 3/4")  
Large: 20 x 17 x 2,5 cm (8 x 6 1/2 x 1")



### Back widening

Create a wider back and distribute pressure from the back profile. Sold as a kit complete with additionally two tension adjustable straps.

Only in combination with seat width 40 cm and 45 cm.



### Back extension 3A

Add 6 cm on 3A- back support profile. With the back extension mounted on to the 3A back support profile a back height of up to 68 cm is achieved.



### Gas piston clips


Every clip reduces range of tilt by 2° or back angle by 2.5°. Useful when a more active balance position is preferred. No tools are required.

Only in combination with Etac Prio manual tilt.



### Positioning belt

With snap-lock, length adjustable. Two alternative attachment points on the frame. Available in two lengths, 100 cm and 128 cm



### Attachment positioning belt, front

Use for front assembly of 4-point belt. Or for alternative placement of 2-point belt. Sold as pair.



### Chest strap attachment

To be combined with head support.



### Electric tilt

Electrical back and seat recline is available as a configuration option or as an accessory. *Can not be combined with Dynamic posterior pelvic support, comfort lock, locking pin, quick release back, gas piston clips and/or basket.*



### Safety edge

To be fitted together with Prio Electrical tilt equipped with 16" transit wheels.



### User operated self-tilt

Enables the more active user to independently manage the adjustment of the seat angle for rest and activity.

To be applied on the seat tilt function only.



### Push handles

Height adjustable and detachable, grips with anti-slide handle covers. Can be combined with attendant brake.



### Bow handle

Height adjustable, detachable. Can be combined with attendant brake.

*Also available as angle adjustable for seat width 40 cm, 45 cm and 50 cm.*



### Attendant brake

With anti-slide handle cover. 10 cm (4") between the handle grip and brake lever.



### Brake with angled handle

Provides a larger gripping surface. Available in black or red. For 20", 22" and 24" rear wheels.



### Brake with extension handle

The brake lever can be angled forward to facilitate sideways transfers.



### One hand brake

One hand brake for right or left side. It can be combined with attendant brakes.

Fits 24" and 22" wheels.



### Lever locking device

Prevents activation of the recline or tilt. To adjust the position, pull the locking device upwards, and then push it back down to lock.



### Comfort Lock

This function ensures the same preset tilt position. There are four alternative settings (0°, 4°, 8° or 12°). Choose one position that is frequently used, for activity or rest. Not in combination with Etac Prio electrical tilt.



### Locking pin quick release back

To facilitate the folding of the back support the original pin can be replaced by quick release solution.

*Not in combination with Etac Prio electrical tilt.*



### Back angle control

Limits the back angle forward on the Prio for an active seat position.



### One arm drive

Fits the left or right sides of the chair. Solid innertube. Vinyl coated handrim.

*Can not be combined with attendant brake, one hand brake, widening kit or basket*



### Additional seat heights for 16" rear wheels

Distance that enables two additional seat heights for 16" rear wheels: 40 cm (15<sup>3</sup>/<sub>4</sub>") and 42.5 cm (16<sup>3</sup>/<sub>4</sub>").



### Rear wheels

20", 22" or 24". Solid and pneumatic tyres. Quick release hub, aluminum, plastic or stainless handrim. 16" transit wheels. Standard wheel for Prio adds only 19 cm to the seat width. With adjustable handrim it adds 21 cm to the seat width.



### Telescopic anti-tips

Foldable, height, adjustable in depth and angle. Fits all seat heights.



### Castor wheels

6", 7" or 8".



### IV-stand

Including all attachments.

*Not in combination with tray.*



### Widening castor housing

Moves the castor 3 cm (1¼") outwards. The chair's total width is not changed.



### Basket for oxygene tank

Can also be used as a carrier basket. The basket can not be combined with Etac Prio electrical tilt.



### Spoke guards

Available in grey or yellow print.



### Tool kit

Contains all necessary tools.  
Allen key 3, 4, 5 and 6 mm  
Ring spanner 8, 10, 13 and 17 mm  
Socket spanner 19/24 mm

## Our Add-On foot plate system encourage movement in your feet

By moving your feet, even if it's small movements, you promote blood circulation and you'll get activation from your hips down, which helps to avoid agitation, pain and discomfort by being too still, while sitting in a wheelchair. Small movements can especially be beneficial to an individual who has always previously been walking.



The Add-Ons with Extended foot supports and Foot box are mounted on to existing foot supports to create a greater surface. The Distal thigh support helps to align the femur, position the knees in a more neutral position and to control abduction. The Add-On foot supports can easily be swung away for easy transfer in and out of the wheelchair. The foot box is fixed, provides padding to decrease pressure and is a gentle aid to help position the feet more comfortably. *The Add-Ons are available as accessories, not in configuration.*



January 4, 2022

Founded in 1973, Etac is a world-leading developer and provider of ergonomic assistive devices and patient handling equipment. Our ambition is to provide solutions that optimise quality of life for the individual, family members and caregivers.

For the latest news and continuously updated product information – visit: [www.etac.com](http://www.etac.com)

Etac AB  
Kista Science Tower  
SE-164 51 Kista, Sweden  
Tel +46 371 - 58 73 30  
[info@etac.se](mailto:info@etac.se) [www.etac.com](http://www.etac.com)

 **etac**<sup>®</sup>  
Creating Possibilities