

F35 - Frequently Asked Questions

- **Does the F35 have a sufficient range for getting around every day?**

With 15km range and powerful 250W motors the F35 is perfect for everyday commuting to the workplace and city breaks.

- **I have shoulder pain, would I still be able to manage the F35?**

The F35 can be adjusted on the subframe to be brought closer or further away to prevent shoulder pain or fatigue.

- **How easy is it to attach and remove from the wheelchair?**

With an easy two movement docking process and quick release the F35 and subframe are light and easy to attach and remove from the chair. The subframe can be left clamped on to the chair when the user is not transferring out of the chair.

- **Can I load and unload this from my car myself?**

With a 7.5kg transport weight and folding handles the F35 is a lightweight and compact package ideal for independent car transfers.

- **Can the F35 handle everyday terrain e.g. going into shops, tight pavements, paving stones?**

The F35 has a solid and durable tyre to tackle everyday terrain, coupled with 180° rotation in the tube to allow the user to turn on the spot indoors and in tight urban environments.

- **Do the brackets need to stay on my wheelchair?**

The F35 clamps can be removed from the chair frame however to ensure correct positioning and easy docking we would recommend keeping the clamps on the chair after they have been fitted correctly with minimal weight added to the chair.

